

TNFA MAB (REFERENCE MAB)

Cat.#: 28064

Product Name: Anti-TNF α Reference Monoclonal Antibody

Synonyms: DIF;NF- α ;NFA;NFSF2;achexin;achectin;NF α

Description: Anti-TNF α Reference Monoclonal Antibody

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Applications: ELISA; Flow Cyt

Recommended Dilutions: ELISA 1:5000-10000; Flow Cyt 1:100

Host Species: Homo sapiens

Isotype: IgG1

Purification: Purified from cell culture supernatant by affinity chromatography

Species Reactivity: Human TNFA

Constituents: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

Storage & Shipping: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. ELISA plate pre-coated by $1 \mu\text{g}/\text{ml}$ ($100 \mu\text{l}/\text{well}$) Human TNF α , His tagged protein 11199 can bind Anti-TNF α Neutralizing antibody (28064) in a linear range of 0.3-500 ng/ml.




Figure 2. Flow cytometry analysis with Anti-TNF α mAb $1 \mu\text{g}/\text{ml}$ (Adalimumab) on Expi293 cells transfected with human TNF α (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).