

CD33 MAB (REFERENCE MAB)

目录: 28024

产品名称: CD33 Reference Monoclonal Antibody

基因符号: CD33;IGLEC3;p67

描述: CD33 Reference Monoclonal Antibody

背景: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

经过测试的应用: Flow Cyt

推荐稀释度: Flow Cyt 1:100

种属反应性: Humanized

亚型: IgG4

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: Human CD33

成分: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form.

After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. HEK293 cell line transfected with irrelevant protein (B) and human CD33 (A) were surface stained with CD33 neutralizing antibody $1\mu\text{g/ml}$ (gemtuzumab) followed by Alexa 488-conjugated human IgG secondary antibody.

