

CD20 MAB (RITUXIMAB)

目录: 28034

产品名称: CD20 (Rituximab) Monoclonal Antibody

基因符号: B1;Bp35;CD20;CVID5;LEU-16;MS4A2;S7

描述: CD20 (rituximab) Monoclonal Antibody

背景: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

经过测试的应用: Flow Cyt

推荐稀释度: Flow Cyt 1:100

种属反应性: Chimeric

亚型: IgG1

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: Human CD20

成分: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. Flow cytometry analysis with 2 $\mu\text{g}/\text{mL}$ CD20 (rituximab biosimilar) mAb (28034) on Expi293 cells transfected with Human CD20 protein (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram). antibody.

