

SAR1A PAB

Sar1a Polyclonal Antibody

货号: 21090

基因符号: SAR1A

描述: Sar1a Rabbit Polyclonal Antibody

背景: Small GTPases are a super-family of cellular signaling regulators. Sar1 is a member of these small GTPases. Sar1 functions as a molecular switch to control protein-protein and protein-lipid interactions that direct vesicle budding from the ER.

免疫原: A synthetic peptide from the internal region of human Sar1a protein

经过测试的应用: ELISA, WB, IHC

推荐稀释度: ELISA 1:100-1:5000

WB 1:100-1:2000

IHC 1:100-1:1000

浓度: 0.2 mg/ml

种属反应性: Rabbit

形式: Liquid

Z43克隆性: Polyclonal

纯化: Purified from serum


Preservative: No

成分: PBS (without Mg^{2+} and Ca^{2+}), pH 7.4, 150 mM NaCl, 50% glycerol

种属反应性: Sar1a antibody recognizes Sar1a proteins of vertebrates.

储存条件: Store at $-20^{\circ}C$. Avoid repeated freezing and thawing

Western blot:



WB: anti-Sar1a pAb

Western blot analysis of Sar1a proteins.

MEF whole-cell lysates were used for Western blot with anti Sar1a rabbit polyclonal antibody (Cat# 21090).