

## **SARS-COV-2 RBD (DM54) RABBIT MAB**

目录: 28353

产品名称: SARS-CoV-2 RBD(DM54) Rabbit Monoclonal Antibody

基因符号: SARS-CoV-2 RBD

描述: SARS-CoV-2 RBD antibody(DM54) Rabbit Monoclonal Antibody

背景: SARS-CoV-2 (Severe Acute Respiratory Syndrome Coronavirus 2) also known as Covid19 (2019 Novel Coronavirus) is a virus that causes illnesses ranging from the common cold to severe diseases. The spike protein is a type I transmembrane protein containing two subunits; S1 and S2. S1 mainly contains a receptor binding domain (RBD); which accounts for recognizing the cell surface receptor; ACE2. S2 contains basic elements needed for the membrane fusion. Recent publications indicate that S1-RBD domain can induce virus neutralizing-antibody and T cell response.

经过测试的应用: ELISA; Flow Cyt

推荐稀释度: ELISA 1:5000-10000; Flow Cyt 1:100

种属反应性: Rabbit

亚型: Rabbit IgG

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: SARS-CoV-2 S protein RBD

成分: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form.

After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. Elisa plate pre-coated by 2  $\mu\text{g/ml}$  (100 $\mu\text{l/well}$ ) SARS-CoV-2 RBD protein can bind Rabbit SARS-CoV-2 RBD monoclonal antibody (clone:DM54) in a linear range of 0.1-100 ng/ml.