

IL6 MAB (SILTUXIMAB)

目录: 28016

产品名称: IL6 (Siltuximab) Monoclonal Antibody

基因符号: IL6;Interleukin-6;BSF2;HSF;IFNB2

描述: IL6 (siltuximab) Monoclonal Antibody

背景: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

经过测试的应用: ELISA; Flow Cyt

推荐稀释度 : ELISA 1:5000-10000; Flow Cyt 1:100

种属反应性: Chimeric

亚型: IgG1

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: Human IL6

成分: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form.

After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).

Anti-IL6 (siltuximab biosimilar) mAb ELISA

0.2 µg of Human IL6, mFc-His Tagged protein per well




Figure 1. ELISA plate pre-coated by 2 µg/ml (100 µl/well) Human IL6, mFc-His tagged protein 11160 can bind IL6 Neutralizing antibody in a linear range of 7.81-386.8 ng/ml.