

SLC7A5 PAB

SLC7A5 Polyclonal Antibody

Cat. #: 21118

Gene Symbol: SLC7A5

Description: Anti-SLC7A5 Rabbit Polyclonal Antibody

Background: Solute carrier family 7 member 5 (also known as large neutral amino acids transporter small subunit 1) is a protein encoded by the SLC7A5 gene in human. SLC7A5 plays a role in neuronal cell proliferation, is involved in the uptake of methylmercury (MeHg) and plays a role in metal ion homeostasis and toxicity.

Immunogen: A synthetic peptide from the internal region of SLC7A5 of human origin.

Applications: ELISA, WB, IHC, IP

Recommended Dilutions:

ELISA: 1:100–1:5000

WB: 1:100–1:2000

IHC: 1:50–1:100

Concentration: 0.25 mg/ml

Host Species: Rabbit

Format: Liquid

Clonality: Polyclonal

Isotype: IgG

Purity: Purified from serum


Preservative: No

Constituents: PBS (without Mg^{2+} and Ca^{2+}), pH 7.4, 150 mM NaCl, 50% glycerol

Species Reactivity: Anti-SLC7A5 antibody recognizes SLC7A5 of vertebrates.

Storage Conditions: Store at $-20^{\circ}C$. Avoid repeated freezing and thawing

Immunohistochemistry:



Western blot analysis of Ric-8A proteins.

Whole-cell lysate of HeLa cells was blotted with Anti-Ric-8A monoclonal antibody (Cat. #26071).