

**SMO(A68V)****SMO(A68V)**

货号: 26425

规格: 100 µL

基因符号: SMO; MRD21

描述: SMO(A68V) 小鼠单克隆抗体

背景: Smoothened is a protein that in humans is encoded by the SMO gene. Smoothened is a G protein-coupled receptor that is a component of the hedgehog signaling pathway and is conserved from flies to humans. It is the molecular target of the teratogen cyclopamine. SMO can function as an oncogene. Activating SMO mutations can lead to unregulated activation of the hedgehog pathway and cancer

免疫原: A synthetic peptide from the internal region of SMO which includes the mutation of A68V, human origin.

经过测试的应用: ELISA, WB, IHC

推荐稀释度:

ELISA: 1:1000-1:5000

WB: 1:500-1:1000

IHC: 1:50-1:100

浓度: 1 mg/ml

种属反应性: Mouse

形式: Liquid

克隆性: Monoclonal

亚型: IgG


纯化: Purified from ascites

**Preservative:** No

成分: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150 mM NaCl, 50% glycerol


种属反应性: Recognizes A68V mutant, but not wild type SMO of vertebrates.

储存条件: Store at -20°C. Avoid repeated freezing and thawing



## Western blot: WB: anti-SMO (A68V) mAb

Western blot analysis of recombinant SMO(A68V) and wild type proteins. Purified His-tagged SMO(A68V) protein (lane 1) and corresponding wild type protein (lane 2) were blotted with SMO(A68V) monoclonal antibody (货号26425).



## Immunofluorescence:

Immunofluorescence of cells expressing SMO proteins with SMO(A68V) antibody. HEK293T cells were transfected with pCDNA3-GFP-SMO (WT) plasmid (left column) or pCDNA3-GFP-SMO(A68V) plasmid (right column), then fixed and stained with SMO(A68V) monoclonal antibody (货号26425).