

## **SSTR2 (DMC428) IGG1 CHIMERIC MAB**

目录: 28222

产品名称: SSTR2(DMC428) IgG1 Chimeric Monoclonal Antibody

基因符号: SRIF-1; SS2R

描述: SSTR2 antibody(DMC428) IgG1 Chimeric Monoclonal Antibody

背景: Somatostatin acts at many sites to inhibit the release of many hormones and other secretory proteins. The biologic effects of somatostatin are probably mediated by a family of G protein-coupled receptors that are expressed in a tissue-specific manner. SSTR2 is a member of the superfamily of receptors having seven transmembrane segments and is expressed in highest levels in cerebrum and kidney.

经过测试的应用: Flow Cyt

推荐稀释度: Flow Cyt 1:100

种属反应性: Rabbit

亚型: Rabbit:Human Fc chimeric IgG1

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: Human SSTR2

成分: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. Flow cytometry analysis with SSTR2 (DMC428) on Expi293 cells transfected with human SSTR2 (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).