

## STK11 MAB

产品名称: STK11 Antibody

货号: 26158

规格: 100  $\mu$ L

基因符号: STK11, LKB1, PJS, hLKB1

描述: STK11 小鼠单抗

背景: Serine/threonine kinase 11 (STK11) also known as liver kinase B1 (LKB1) or renal carcinoma antigen NY-REN-19 is a kinase that controls the activity of AMP-activated protein kinase (AMPK) family members, thereby playing a role in various processes such as cell metabolism, cell polarity, apoptosis and DNA damage response, and functions as a tumour suppressor. Recent studies have uncovered a large number of somatic mutations of the STK11 gene that are present in lung, cervical, breast, intestinal, testicular, pancreatic and skin cancer.

免疫原: A synthetic peptide from the internal region of STK, which includes the mutation of P281L, human origin.

经过测试的应用: ELISA, WB, IHC

推荐稀释比:

ELISA: 1:1000-1:5000

WB: 1:100-1:1000

IHC: 1:50-1:200

**Concentration:** 0.5 mg/ml

种属反应性: Mouse

形式: Liquid

克隆性: Monoclonal


亚型: IgG

纯化: Purified from ascites

成分: PBS (without  $Mg^{2+}$  and  $Ca^{2+}$ ), pH 7.4, 150 mM NaCl, 50% glycerol

种属反应性: Recognizes P281L mutated, but not wild type STK11 of vertebrates.

储存条件: Store at  $-20^{\circ}C$ . Avoid repeated freezing and thawing



### Immunofluorescence:

Immunofluorescence of cells expressing STK11 proteins with Anti-STK11 (P281L) antibody. HEK293T cells were transfected with pCDNA3-GFP-STK11 (WT) plasmid (left column) or pCDNA3-GFP-STK11 (P281L) plasmid (right column), then fixed and stained with Anti-STK11 (P281L) monoclonal antibody (货号:26158).