

GPRC5D MAB (TALQUETAMAB)

Cat.#: 28075

Product Name: Anti-GPRC5D(Talquetamab) Monoclonal Antibody

Synonyms: GPRC5D

Description: Anti-GPRC5D(Talquetamab) Monoclonal Antibody

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Applications: Flow Cyt

Recommended Dilutions: Flow Cyt 1:10

Host Species: Humanized

Isotype: IgG4

Purification: Purified from cell culture supernatant by affinity chromatography

Species Reactivity: Human GPRC5D

Constituents: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

Storage & Shipping: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. Flow cytometry analysis with Anti-GPRC5D mAb 25 $\mu\text{g}/\text{ml}$ (Talquetamab) on Expi293 cells transfected with human GPRC5D (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram) .



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
