

## NOTCH3 MAB (TAREXTUMAB)

**Cat.#:** 28099

**Product Name:** Anti-NOTCH3(Tarextumab) Monoclonal Antibody

**Synonyms:** CADASIL;CADASIL1;CASIL;IMF2;LMNS

**Description:** Anti-NOTCH3(tarextumab ) Monoclonal Antibody

**Background:** Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

**Applications:** Flow Cyt

**Recommended Dilutions:** Flow Cyt 1:100

**Host Species:** Homo sapiens

**Isotype:** IgG2

**Purification:** Purified from cell culture supernatant by affinity chromatography

**Species Reactivity:** Human NOTCH3

**Constituents:** Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

**Storage & Shipping:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. Flow cytometry analysis of MCF-7 cells stained with 1 µg/mL Anti-NOTCH3 (tarextumab biosimilar) mAb (28099) (Blue histogram) and unstained MCF-7 cells as a negative control (Red histogram).