

NOTCH3 MAB (TAREXTUMAB)

目录: 28099

产品名称: NOTCH3(Tarextumab) Monoclonal Antibody

基因符号: CADASIL;CADASIL1;CASIL;IMF2;LMNS

描述: NOTCH3(tarextumab) Monoclonal Antibody

背景: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

经过测试的应用: Flow Cyt

推荐稀释度: Flow Cyt 1:100

种属反应性: Homo sapiens

亚型: IgG2

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: Human NOTCH3

成分: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form.

After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. Flow cytometry analysis of MCF-7 cells stained with 1 µg/mL NOTCH3 (tarextumab biosimilar) mAb (28099) (Blue histogram) and unstained MCF-7 cells as a negative control (Red histogram).