

NKG2D MAB (TESNATILIMAB)

Cat.#: 28048

Product Name: Anti-NKG2D (Tesnatilimab) Monoclonal Antibody

Synonyms: NKG2D;D314;LRK1;K cell receptor D

Description: Anti-NKG2D (tesnatilimab) Monoclonal Antibody

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Applications: ELISA; Flow Cyt

Recommended Dilutions: ELISA 1:5000-10000; Flow Cyt 1:100

Host Species: Homo sapiens

Isotype: IgG4

Purification: Purified from cell culture supernatant by affinity chromatography

Species Reactivity: Human NKG2D

Constituents: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

Storage & Shipping: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).

Anti-NKG2D (tesnatilimab biosimilar) mAb ELISA

0.2 µg of Human NKG2D, mFc Tagged protein per well




Figure 1. ELISA plate pre-coated by 2 µg/ml (100 µl/well) Human NKG2D, mFc tagged protein 11200 can bind Anti-NKG2D Neutralizing antibody (28048) in a linear range of 0.64-400 ng/ml.




Figure 2. Flow cytometry analysis with 15 µg/mL Anti-NKG2D (tesnatilimab) mAb (28048) on Expi293 cells transfected with Human NKG2D protein and Human DAP10 protein (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).