

## TFF3 MAB

### TFF3 Monoclonal Antibody

货号: 26047

基因符号: TFF3

描述: TFF3 小鼠单克隆抗体

背景: TFF3 (trefoil factor 3) is a member of the trefoil family proteins which are characterized by having at least one copy of the trefoil motif, a 40-amino acid domain that contains three conserved disulfides. It is involved in the maintenance and repair of the intestinal mucosa, and promotes the mobility of epithelial cells in healing processes. TFF3 is part of a panel of urinary biomarkers to monitor reversibility of renal injury.

免疫原: A synthetic peptide from the internal region of TFF3, human origin.

经过测试的应用: ELISA, WB, IHC, IP

推荐稀释度:

ELISA: 1:1000-1:5000

WB: 1:500-1:1000

IHC: 1:50-1:100

浓度: 1 mg/ml

种属反应性: Mouse

形式: Liquid

克隆性: Monoclonal

亚型: IgG

纯化: Purified from ascites


**Preservative:** No

成分: PBS (without  $Mg^{2+}$  and  $Ca^{2+}$ ), pH 7.4, 150 mM NaCl, 50% glycerol

种属反应性: TFF3 monoclonal antibody recognizes TFF3 of vertebrates.

储存条件: Store at  $-20^{\circ}C$ . Avoid repeated freezing and thawing

## Western blot:



Western blot analysis of recombinant TFF3 proteins. Purified GST-tagged TFF3 proteins (货号10159) were blotted with anti TFF3 monoclonal antibody (Cat # 26047).