

SSTR2 MAB (TIDUTAMAB)

Cat.#: 28122

Product Name: Anti-SSTR2(Tidutamab) Monoclonal Antibody

Synonyms: SS-2-R;SS2-R;SS2R;SST2;SRIF-1

Description: Anti-SSTR2(tidutamab) Monoclonal Antibody

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Applications: ELISA; Flow Cyt

Recommended Dilutions: ELISA 1:5000-10000; Flow Cyt 1:100

Host Species: Chimeric Humanized

Isotype: IgG1

Purification: Purified from cell culture supernatant by affinity chromatography

Species Reactivity: Human SSTR2

Constituents: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

Storage & Shipping: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. ELISA plate pre-coated by 2 µg/mL (100 µl/well) Human SSTR2 Protein, hFc Tag (11591) can bind Anti-SSTR2 Neutralizing antibody (28122) in a linear range of 2.93-187.50 ng/mL. In order to specifically detect 28122, mouse anti-human Fab-specific antibody was used as detection antibody.




Figure 2. Flow cytometry analysis with 1 µg/mL Anti-SSTR2 (tidutamab biosimilar) mAb (28122) on Expi293 cells transfected with Human SSTR2 protein (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).