

HUMAN TLR9 FULL LENGTH PROTEIN

目录: 11127

产品名称: Human TLR9 Full Length Protein

规格: 10 µg, 50 µg and 100 µg

基因符号: CD289

Target: TLR9

UNIPROT ID: Q9NR96

描述: Human TLR9 full length protein-synthetic nanodisc

背景: The protein is a member of the Toll-like receptor (TLR) family, which plays a fundamental role in pathogen recognition and activation of innate immunity. TLRs are highly conserved from *Drosophila* to humans and share structural and functional similarities. They recognize pathogen-associated molecular patterns (PAMPs) that are expressed on infectious agents, and mediate the production of cytokines necessary for the development of effective immunity. Studies in mice and human indicate that this receptor mediates cellular response to unmethylated CpG dinucleotides in bacterial DNA to mount an innate immune response.

Species/Host: HEK293

Molecular Weight: The human full length TLR9 protein has a MW of 115.9 kDa

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure1. Elisa plates were pre-coated with Flag Tag TLR9-Nanodisc (0.2 µg/per well). Serial diluted Flag monoclonal antibody solutions were added, washed, and incubated with secondary antibody before Elisa reading. From above data, the EC50 for Flag monoclonal antibody binding with TLR9-Nanodisc is 2.467ng/ml.




Figure2. Human TLR9-Nanodisc, Flag Tag on SDS-PAGE