

TP53(R282W) PAB

TP53(R282W) Polyclonal Antibody

货号: 21165

基因符号: TP53

描述: TP53(R282W) Rabbit Polyclonal Antibody

背景: p53 protein acts as a tumor suppressor in many tumor types and induces growth arrest or apoptosis depending on the physiological circumstances and cell type. Defects in TP53 are found in many types of cancers.

免疫原: A synthetic peptide near the C-terminus of human TP53 which includes the mutation of R282W.

经过测试的应用: ELISA, WB, IHC, IF

推荐稀释度:

ELISA: 1:1000-1:5000

WB: 1:1000-1:5000

IHC: 1:100-1:1000

浓度: 0.1 mg/ml

种属反应性: Rabbit

形式: Liquid

克隆性: Polyclonal

亚型:

纯化:: Purified from serum

Preservative: No

成分: PBS (without Mg^{2+} and Ca^{2+}), pH 7.4, 150 mM NaCl, 50% glycerol


种属反应性: TP53(R282W) antibody recognizes TP53(R282W) of vertebrates.

储存条件:: Store at $-20^{\circ}C$. Avoid repeated freezing and thawing.

Western blot:

TP53 protein

WT R282W



WB: Anti-TP53(R282W) pAb

Western blot analysis of recombinant TP53(R282W) and wild type proteins. Purified His-tagged TP53(R282W) protein(目录 10204) and wild type protein (货号 10202) were blotted with TP53(R282W) rabbit antibody (货号 21165).