

## OSMR MAB (VIXARELIMAB)

目录: 28140

产品名称: OSMR(Vixarelimab) Monoclonal Antibody

基因符号: IL-31 receptor subunit beta;IL-31R subunit beta;IL-31R-beta;IL-31RB

描述: OSMR(vixarelimab) Monoclonal Antibody

背景: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

经过测试的应用: Flow Cyt

推荐稀释度: Flow Cyt 1:100

种属反应性: Homo sapiens

亚型: IgG4-lambda2

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: Human OSMR

成分: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form.

After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. Flow cytometry analysis with 15  $\mu\text{g}/\text{mL}$  OSMR (vixarelimab biosimilar) mAb (28140) on Expi293 cells transfected with Human OSMR protein (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).