

## Product Description

### Rab3A

**Catalog Number:** 21041

**Gene Symbol:** RAB3A

**Description:** Anti-Rab3A Rabbit Polyclonal Antibody

**Background:** Rab3A is a member of the Rab family GTPases and is specifically targeted to synaptic vesicles.

**Immunogen:** A synthetic peptide near the C-terminus of human Rab3A.

**Tested Applications:** ELISA, WB, IHC, IP

**Recommended Dilutions:**

ELISA 1:100-1:5000

WB 1:100-1:2000

IHC 1:50-100

**Concentration:** 0.2 mg/ml

**Host:** Rabbit

**Clonality:** Polyclonal

**Purity:** Purified from serum

**Format:** Liquid


**Preservative:** No

**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH7.4, 150 mM NaCl, 50% glycerol

**Species Reactivity:** Rab3A antibody recognizes Rab3A of vertebrates.

**Storage Conditions:** Store at -20°C. Avoid freeze / thaw cycles

**Western blot:**



**WB: anti-Rab3A pAb**

**Western blot analysis of anti-Rab3A antibody.**

Purified GST-tagged Rab3A protein was used for Western blot with anti-Rab3A polyclonal antibody (Cat. # 21041).

FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC APPLICATIONS