

Rab25

Catalog Number: 21056

Gene Symbol: RAB25

Description: Anti-Rab25 Rabbit Polyclonal Antibody

Background: Ras-related protein Rab25 is involved in recycling vesicles. It has been associated with the progression and aggressiveness of certain tumors.

Immunogen: A synthetic peptide from the internal region of human Rab25 protein

Tested Applications: ELISA, WB, IHC

Recommended Dilutions:

ELISA 1:100-1:5000

WB 1:100-1:1000

IHC 1:100-1:1000

Concentration: 0.2 mg/ml

Host: Rabbit

Clonality: Polyclonal

Purity: Purified from serum

Format: Liquid


Preservative: No

Constituents: PBS (without Mg^{2+} and Ca^{2+}), pH7.4, 150 mM NaCl, 50% glycerol

Species Reactivity: Rab25 antibody recognizes Rab25 proteins from vertebrates.

Storage Conditions: Store at $-20^{\circ}C$. Avoid freeze / thaw cycles

Western blot:



WB: anti-Rab25 pAb

HEK-293T whole-cell lysate was used for Western blot with anti-Rab25 polyclonal antibody (Cat. # 21056).

FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC APPLICATIONS