

Rab12

Catalog Number: 21076

Gene Symbol: RAB12

Description: Anti-Rab12 Rabbit Polyclonal Antibody

Background: Rab12 protein is a member of the Rab subfamily of Ras-related small GTPases. Rab12 might regulate vesicular trafficking.

Immunogen: A synthetic peptide from the internal region of Rab12

Tested Applications: ELISA, WB, IHC, IP

Recommended Dilutions:

ELISA	1:100-1:5000
WB	1:100-1:2000
IHC	1:10-100

Concentration: 0.2 mg/ml

Host: Rabbit

Clonality: Polyclonal

Purity: Purified from serum

Format: Liquid


Preservative: No

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH7.4, 150 mM NaCl, 50% glycerol

Species Reactivity: Rab12 antibody recognizes Rab12 from vertebrates.

Storage Conditions: Store at -20°C. Avoid freeze / thaw cycles

Western blot:



HEK 293T cell lysate was used for Western blot with anti-Rab12 polyclonal antibody (Cat. # 21076).

FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC APPLICATIONS