

Product Description

Pioneering GTPase and Oncogene Product Development since 2010

B-RAF

Catalog Number: 26129

Gene Symbol: BRAF, BRAF1, RAFB1

Description: Anti-B-Raf Mouse Monoclonal Antibody

Background: B-Raf belongs to the TKL Ser/Thr protein kinase family, RAF subfamily. It is involved in the transduction of mitogenic signals from the cell membrane to the nucleus and may play a role in the postsynaptic responses of hippocampal neuron. B-Raf is located to nucleusand cytoplasm and expressed in brain and testis.

Immunogen: A synthetic peptide from the internal

region of B-Raf, human origin.

Tested Applications: ELISA, WB, IF, IHC

Recommended Dilutions:

ELISA: 1:1000-1:5000 WB: 1:100-1:1000 IHC: 1:50-1:100

Concentration: 1 mg/ml

Host: Mouse

Clonality: Monoclonal

Isotype: IgG

Purity: Purified from ascites

Format: Liquid Preservative: No

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH7.4, 150

mM NaCl, 50% glycerol


Species Reactivity: Anti-B-Raf antibodyrecognizes B-

Raf of vertebrates.

Storage Conditions: Store at -20°C. Avoid repeated

freezing and thawing

Western blot:



Western blot analysis of recombinant B-RAF protein.

Purified His-tagged B-Raf protein (Cat. # 10161) was blotted with anti-B-Raf monoclonal antibody (Cat. #

26129).

info@neweastbio.com