

PIK3CA (Q546R)

Catalog Number: 26214

Gene Symbol: PIK3CA, CLOVE, CWS5, MCAP, MCM, MCMTC, PI3K, p110-alpha

Description: Anti-PIK3CA (Q546R) Mouse Monoclonal Antibody

Background: Phosphatidyl 3-kinases (PI3K) are a family of lipid kinases involved in many cellular processes, including cell growth, proliferation, differentiation, motility, and survival. It is composed of an 85 kDa regulatory subunit and a 110 kDa catalytic subunit. Mutants of PIK3CA has been found in the pathogenesis of several cancers, including colon cancer, gliomas, gastric cancer, breast cancer, endometrial cancer, and lung cancer.

Immunogen: A synthetic peptide from the internal region of PIK3CA, which includes the mutation of Q546R, human origin.

Tested Applications: ELISA, WB

Recommended Dilutions:

ELISA: 1:1000-1:5000

WB: 1:100-1:1000

Concentration: 1 mg/ml

Host: Mouse

Clonality: Monoclonal

Isotype: IgG

Purity: Purified from ascites

Format: Liquid


Preservative: No

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH7.4, 150 mM NaCl, 50% glycerol

Species Reactivity: Recognizes Q546R mutant, but not wild-type PIK3CA of vertebrates.

Storage Conditions: Store at -20°C. Avoid repeated freezing and thawing

Western blot:



Western blot analysis of recombinant PIK3CA (Q546R) and wildtype proteins.

Purified His-tagged PIK3CA (Q546R) protein (lane 2) and corresponding wildtype protein (lane 1) were blotted with anti-PIK3CA (Q546R) monoclonal antibody (Cat. #26214).

FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC APPLICATIONS