

## PTEN(R130G)

**Catalog Number:** 26312

**Gene Symbol:** PTEN, 10q23del, BZS, DEC, GLM2, MHAM, MMAC1, PTEN1, TEP1, MMAC1

**Description:** Anti-PTEN(R130G) Mouse Monoclonal Antibody

**Background:** Phosphatase and tensin homolog (PTEN) was identified as a tumor suppressor that is mutated in a large number of cancers at high frequency. It is one of the most commonly lost tumor suppressors in human cancer. During tumor development, mutations and deletions of PTEN occur that inactivate its enzymatic activity leading to increased cell proliferation and reduced cell death. Frequent genetic inactivation of PTEN occurs in glioblastoma, endometrial cancer, and prostate cancer; and reduced expression is found in many other tumor types such as lung and breast cancer.

**Immunogen:** A synthetic peptide from the internal region of PTEN, which includes the mutation of R130G, human origin.

**Tested Applications:** ELISA, WB

**Recommended Dilutions:**

ELISA: 1:1000-1:2000

WB: 1:500-1:1000

**Concentration:** 1.0 mg/ml

**Host:** Mouse

**Clonality:** Monoclonal

**Isotype:** IgG

**Purity:** Purified from ascites

**Format:** Liquid


**Preservative:** No

**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH7.4, 150 mM NaCl, 50% glycerol

**Species Reactivity:** Recognizes PTEN(R130G), but not wild type PTEN protein from vertebrates.

**Storage Conditions:** Store at -20°C. Avoid freeze / thaw cycles.

### Western blot:



### Western blot analysis of recombinant PTEN (R130G) and wild type proteins.

Purified His-tagged FGFR1 (R130G) protein (lane 2) and corresponding wild type protein (lane 1) were blotted with anti-PTEN (R130G) monoclonal antibody (Cat. # 26312).

FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC APPLICATIONS