

## Ras(G13R)

**Catalog Number:** 26474

**Gene Symbol:** H-Ras; K-Ras; N-Ras

**Description:** Anti-RAS(G13R) Mouse Monoclonal Antibody

**Background:** Ras protein is a membrane-associated GTPase signaling protein that regulates proliferation, differentiation, and cell survival. Ras G13R mutation result in decreased GTPase activity and constitutive signaling and is correlated with the development of pancreatic intraepithelial neoplasia and lung lesions.

**Immunogen:** A synthetic peptide from the internal region of Ras which includes the mutation of G13R, human origin.

**Tested Applications:** ELISA, WB, IF, IHC

**Recommended Dilutions:**

ELISA: 1:1000-1:5000

WB: 1:500-1:1000

IF: 1:50-1:100

IHC: 1:20-1:100

**Concentration:** 1 mg/ml

**Host:** Mouse

**Clonality:** Monoclonal

**Purity:** Purified from ascites

**Format:** Liquid


**Preservative:** No

**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH7.4, 150 mM NaCl, 50% glycerol

**Species Reactivity:** Recognizes G13R mutant, but not wild-type RAS of vertebrates.

**Storage Conditions:** Store at -20°C. Avoid freeze / thaw cycles.

**Western blot:**




**Western blot analysis of recombinant HRas (G13R) and wild type proteins.**

Purified His-tagged Ras (G13R) protein (lane 1) and corresponding wild type protein (lane 2) were blotted with anti-Ras (G13R) monoclonal antibody (Cat. # 26474).

FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC APPLICATIONS

## Immunofluorescence:



## Immunofluorescence of cells expressing Ras proteins with anti-Ras (G13R) antibody.

HEK293T cells were transfected with pCDNA3-GFP-Ras (G13R) plasmid, pCDNA3-GFP-Ras (WT) plasmid, pCDNA3-GFP-Ras (G13S) plasmid, pCDNA3-GFP-Ras (G13C) plasmid, pCDNA3-GFP-Ras (G13V) plasmid, pCDNA3-GFP-Ras (G13A) plasmid, or pCDNA3-GFP-Ras (G13D) plasmid, then fixed and stained with anti-Ras (G13R) monoclonal antibody (Cat. # 26474).

FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC APPLICATIONS