

EGFR(T854A)

Catalog Number: 26481

Gene Symbol: EGFR/HER1/ERBB1

Description: Anti-EGFR(T854A) Mouse Monoclonal Antibody

Background: EGFR is a receptor tyrosine kinase binding ligands of the EGF family and activating several signaling cascades to convert extracellular cues into appropriate cellular responses. It activates at least 4 major downstream signaling cascades including the RAS-RAF-MEK-ERK, PI3 kinase-AKT, PLCgamma-PKC and STATs modules. Defects in EGFR are associated with lung cancer (LNCR).

Immunogen: A synthetic peptide from the internal region of EGFR, which includes the mutation of T854A, human origin.

Tested Applications: ELISA, WB, IF, IHC

Recommended Dilutions:

ELISA:	1:1000-1:5000
WB:	1:1000-1:5000
IHC:	1:50-1:200

Concentration: 1 mg/ml

Host: Mouse

Clonality: Monoclonal

Purity: Purified from ascites

Format: Liquid


Preservative: No

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH7.4, 150 mM NaCl, 50% glycerol

Species Reactivity: Recognizes EGFR (T854A), but not wild type EGFR from vertebrates.

Storage Conditions: Store at -20°C. Avoid repeated freezing and thawing

Western blot:




Western blot analysis of recombinant EGFR (T854A) and wild type proteins.

Purified His-tagged EGFR (T854A) protein (lane 1) and corresponding wild type protein (amino acids 702-987, lane 2) were blotted with anti-EGFR (T854A) mouse monoclonal antibody (Cat. # 26481).

FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC APPLICATIONS

Immunofluorescence:



Immunofluorescence of cells expressing EGFR proteins with anti-EGFR (T854A) antibody.

HEK293T cells were transfected with pCDNA3-GFP-EGFR (WT) plasmid (left column) or pCDNA3-GFP-EGFR (T854A) plasmid (right column), then fixed and stained with anti-EGFR (T854A) monoclonal antibody (Cat. # 26481).

FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC APPLICATIONS